

aCart™ Compact

DATASHEET

aCart™ Compact is a modern, mobile workstation that adds value to the Ambu single-use product portfolio. The aCart Compact workstation provides convenient storage and portability without taking up a lot of space.

A modern mobile workstation

The aCart Compact workstation enables you to move fast when you need to respond to an emergency or change clinical environments. It is also ideal for the outpatient clinic or for doctor's offices with limited space.

Robust and portable

The design and choice of materials make the aCart Compact workstation light enough to easily transport, yet robust and stable enough to protect your valuable endoscopic equipment. The hygiene of the aCart Compact workstation is certified by the DGKH (German Society for Hospital Hygiene).



KEY BENEFITS

- Height adjustable IV pole
- Ergonomic design handle/ cable winder
- Monitor/displaying and processing unit holder
- Preparation tray/shelve
- Storage hook for spare endoscopes
- Support column with integrated cable duct
- Storage basket for accessories, cables or airway tools
- Base with four lockable double swivel casters

Ambu

TECHNICAL SPECIFICATIONS

Dimensions [mm (in)]

Width	530 (20.86)
Height without IV pole	1320 (51.97)
Height with IV pole min./max.	1511 (59.49) / 2351 (92.55)
Depth	610 (24.02)
Package Dimensions W/D/H	1180/590/225 (46.46/23.23/8.86)

Weight [kg (lbs)]

Total	23 (50.71)
-------	------------

Load Capacity [kg (lbs)]

Maximum load	35 (77.1)
Shelve	10 (22)
Monitor holder	14 (30.8)
Basket	3 (6.6)
IV Pole	5 (11)
Optional Vario Mount	8 (17.6)

Device Classification

Europe	Class I, CE certified according to DIN EN 60601-1
USA	Class I, 510(k) exempt from premarket notification

Operating environment

Operating temperature	10 - 40°C (50 - 104°F)
Operating relative humidity	30 - 75%
Operating atmospheric pressure	70-106kPa

Storage environment

Storage temperature	-25 - 70°C (-13 - 158°F)
Storage relative humidity	10 - 95%
Storage atmospheric pressure	50-120kPa



ORDERING SPECIFICATIONS

Description

Reference number

aCart Compact	KU.9713.800
aCart Compact Vario Mount	HA.3000.991

Ambu

Distributor:

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com

Manufacturer:

ITD GmbH
Sportplatzstraße 3
84381 Johanniskirchen
Germany