

DATASHEET

The Basic Buddy manikins are lightweight, portable and suitable for an economical and realistic CPR training.

They feature a realistic head tilt and chest rise, palpable landmarks, a sanitary airway system and a lightweight and compact design.

Adult and Child

The Basic Buddy can be used for adult and child training by changing the compression resistance with a foam pad

Sanitary system

The disposable lung/mouth protection bags are changed with a lung insertion tool after each student preventing cross-contamination and the need to sanitize the manikins.

Easy AED training

The closed and anatomically correct torsi allow the placement of defibrillation pads for realistic, simulated AED training.

Note: Use only AED pads designed for foam body manikins.

Basic Buddy®



KEY BENEFITS

- Practice adult and child CPR techniques
- Lightweight, portable and economical and Tough and durable
- Airways open only with correct hyperextension
- Visible chest rise during ventilation
- Rebreathing airway
- Sanitary mouth/lung protection bags
- Easy AED training

ORDERING SPECIFICATIONS

Item no.

Item no.	
LF03694U	Basic Buddy 5-Pack 5 adult/child torsi, 50 lung/mouth protection bags, 5 lung insertion tools, carry case, instruction manual
LF03696U	Basic Buddy Lung Bags, disposable (set of 100)

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com