

## DATASHEET

AmbuMan® Compression iQF is a first aid and training manikin.

And now with the new iQF App, you can provide training with a high level of precision for compression in adults.

### Compressions can be life saving


Compressions help to maintain perfusion and contractile function of the heart at the brain and coronary level.

Therefore, it is important that people are trained in this discipline in order to save lives.

### iQF, Training that Works for Life

Helping health instructors to train. The iQF application includes data-driven information of compression parameters.

Help your students perform perfect compression depth while enjoying games and competing with each other.



Ambu iQF Instructor

### Correct compression and AED Training

The AmbuMan Compression enables students to train on external chest compression which is one of the main focuses of modern BLS algorithms. The compression indicator confirms the correct compression depth. Additionally, the closed shaped torso recreates realistic training of the placement of AED pads.

# AmbuMan® Compression iQF



## KEY BENEFITS

- Robust and reliable cost efficient BLS training manikin with anatomically correct closed Torso
- Visual feedback built into the manikin torso provides instant feedback of compression amplitude
- Easy AED Training
- Comes in a pack of 5 for easy class training

## TRAIN IQF, TRAIN FOR LIFE

Visualization in detail of the compression depth at your fingertips with Instructor and Trainee iQF App. Connected by Bluetooth®.

# Ambu

## SPECIFICATIONS

Weight Torso	2 kg
Weight 5 x Torso including carrying bag	10.5 kg
Length Torso	45 cm
Ma. compression depth	7 cm


## MATERIALS

Basic unit	ABS (acrylonitrile butadiene styrene)
Chest skin	SEBS (styrene-ethylene-butylene-styrene)


## ORDERING SPECIFICATIONS

Item no.	
A268011000	AmbuMan Compression 5 pack iQF


\* Bluetooth is a registered trademark of Bluetooth SIG, inc.



Placement of AED pads



Compression depth indicator



Ambu iQF Trainee

# Ambu

■ Ambu A/S  
Baltorpbakken 13  
2750 Ballerup  
Denmark  
T +45 72 25 20 00  
F +45 72 25 20 50  
ambu.com