

## DATASHEET

---

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope™ Duodeno, you have a new scope with consistent performance every time, and there is no risk of patient cross-contamination. What's more, there's no need for costly reprocessing or repair.



Distal tip showing elevator in use with a sphincterotome



Elevator control lever to raise (up) and lower (down) the elevator

**Ambu® aScope™ Duodeno  
Sterile.  
Single-use.  
Problem solved.**



## KEY BENEFITS

---

Sterile

No risk of patient cross-contamination

Cost transparency

Eliminates the risk and cost of reprocessing

Familiar design and function to that of conventional duodenoscopes

**Ambu**

## SPECIFICATIONS

### Optical System

Direction of View	Retro Side Viewing 6°
Field of View	130°

### Insertion Section

Distal End Outer Diameter	φ 13.7 mm
Insertion Tube Outer Diameter	φ 11.3 mm
Working Length	1240 mm

### Instrument Working Channel

Inner Diameter	φ 4.2 mm (12.6 FR)
----------------	--------------------

### Bending Section

Angulation Range	Up 120° / Down 90° Left 90° / Right 110°
------------------	---

### Illumination

Illumination source	Dual LEDs
---------------------	-----------

### Dimensions


Weight	700 g
--------	-------

## ORDERING SPECIFICATIONS\*

Reference number	Description	Quantity
485001000	aBox Duodeno	1
482001000	aScope Duodeno	Box of 6
482001000/4	aScope Duodeno	Box of 4**

\*Some configurations may not be available worldwide

\*\*The outer box contains four inner boxes, each containing one scope.



# Ambu

Ambu A/S  
 Baltorpbakken 13  
 2750 Ballerup  
 Denmark  
 T +45 72 25 20 00  
 ambu.com

  
 2797

US: Rx Only